ISTRUZIONI D'USO PUTTY SILICONE

CONFEZIONI

cod. 1911903 5000 gr. A + 5000 gr. B
cod. 1911803 1000 gr A + 1000 gr. B

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INSTRUCTIONS FOR USE PUTTY SILICONE

PACKAGES:

cod. 1911903 5000 gr. A + 5000 gr. B
cod. 1911803 1000 gr A + 1000 gr. B

High viscosity polyaddition silicone for professional dental use.

ADVANTAGES:

- Excellent reproduction of details
- Good mechanical properties
- Low linear shrinkage
- No release of toxic substances during and after vulcanization
- High resistance to inorganic chemicals and UV rays
- Good release properties with organic resins and hydrated gypsum

TECHNICAL DATA:

Before catalysis

Appearance: paste, tasteless and odorless
Color: whitish (A) - light green (B)

In the catalysis phase

1 part of Putty A + 1 part of Putty B
Processing time: about 1 min
Stabilization time: ca 2 min
Total setting time: approx. 3 min

After crosslinking

on specimen with a thickness of 6 mm
Hardness after 24 hours: 70 shore
Volumetric mass at 23 ° C: 1.55
Linear shrinkage: <0.2%
Elastic recovery: 99.4%
Compatibility with plaster: excellent

Note: higher temperatures reduce these times, lower temperatures lengthen them.

USE

Mix thoroughly until you get a homogeneous mixture. Attention, it can be inhibited by contact with the following products: heavy mineral salts, amines, sulfur and derivatives, epoxy resin catalysts.

STORAGE

Keep packages sealed at temperatures <23 ° C.